MDSe User Guide

CE

MDSe User Guide

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English

Deutsch

Español

Français

Italiano

Português

SonoSite, Inc.

21919 30th Drive SE

Bothell, WA 98021

USA

T: 1-888-482-9449 or 1-425-951-1200

F: 1-425-951-1201

SonoSite Ltd

Alexander House

40A Wilbury Way

Hitchin

Herts SG4 0AP

UK

T: +44-1462-444800

F: +44-1462-444801

Caution: | United States federal law restricts this device to sale by or on the order of a physician.

"MicroMaxx" is a registered trademark or trademark of SonoSite, Inc.

Non-SonoSite product names may be registered trademarks or trademarks of their respective owners.

SonoSite products may be covered by one or more of the following U.S. patents: 4454884, 4462408, 4469106, 4474184, 4475376, 4515017, 4534357, 4542653, 4543960, 4552607, 4561807, 4566035, 4567895, 4581636, 4591355, 4603702, 4607642, 4644795, 4670339, 4773140, 4817618, 4883059, 4887306, 5016641, 5050610, 5095910, 5099847, 5123415, 5158088, 5197477, 5207225, 5215094, 5226420, 5226422, 5233994, 5255682, 5275167, 5287753, 5305756, 5353354, 5365929, 5381795, 5386830, 5390674, 5402793, 5,423,220, 5438994, 5450851, 5456257, 5471989, 5471990, 5474073, 5476097, 5479930, 5482045, 5482047, 5485842, 5492134, 5517994, 5529070, 5546946, 5555887, 5603323, 5606972, 5617863, 5634465, 5634466, 5636631, 5645066, 5648942, 5669385, 5706819, 5715823, 5718229, 5720291, 5722412, 5752517, 5762067, 5782769, 5800356, 5817024, 5833613, 5846200, 5860924, 5893363, 5916168, 5951478, 6036643, 6102863, 6104126, 6113547, 6117085, 6142946, 6203498 B1, 6371918, 6135961, 6364839, 6383139, 6416475, 6447451, 6471651, 6569101, 6575908, 6604630, 6648826, 6835177, 6817982, 6730035, 6962566, 7169108, D0280762, D0285484, D0286325, D0300241, D0306343, D0328095, D0369307, D0379231, D456509, D461895, D509900. Other patents pending.

P06427-03 04/2007

Copyright 2007 by SonoSite, Inc.

All rights reserved.

MDSe User Guide

Contents

Introduction	1
Safety	2
MDSe Features	4
Connecting Ultrasound System	6
Connecting Printer and DVD Recorder or VCR	6
Connecting Printer Only	11
Connecting External Display	13
Connecting Triple Transducer Connect	16
Cables, Connectors, and Mini-Dock Symbols	17
Troubleshooting	23
Cleaning and Disinfecting the MDSe	25
Cleaning the Peripherals	25
Specifications	26

Introduction

The Mobile Docking System enhanced (MDSe) provides a platform for the MicroMaxx® ultrasound system with connections for accessories and peripherals. This user guide provides information on the MDSe and on connecting accessories and peripherals.

For information on the peripherals, please see the manufacturer's operating instructions shipped with each product. For information on the ultrasound system, transducers, and additional safety information including warnings and cautions, see the ultrasound system user guide.

WARNING: To avoid the risk of patient injury, misdiagnosis, or user injury, read all the warnings in the ultrasound system user guide and supplemental guides.

Conventions

These conventions are used in this user guide:

- A **WARNING** describes precautions necessary to prevent injury or loss of life.
- A **Caution** describes precautions necessary to protect the products.
- Numbered steps in procedures must be performed in a specific order.
- Bulleted lists present information in list format but do not imply a sequence.

Labeling Symbols and Terms

Labeling symbols and terms used on the system are explained in the ultrasound system user guide.

Customer Assistance and Comments

Questions and comments are encouraged. SonoSite is interested in your feedback regarding the system and the user guide. Please call SonoSite at **1-888-482-9449**. For technical support, call **1-877-657-8118**, or if you are outside the USA, call the nearest SonoSite representative. You can also e-mail SonoSite at **comments@sonosite.com**.

Safety

For maximum safety, observe the following warnings and cautions before using the accessories and peripherals with the ultrasound system and transducers.

Electrical Safety

- WARNING:** To avoid the risk of electrical shock and fire hazard, a multiple portable socket outlet (MPSO) or extension cord should not be connected to the auxiliary mains socket on the MDSe.
- To avoid the risk of electrical shock to the patient, do not touch the ungrounded connectors on the back of the mini-dock while simultaneously touching a patient.
- To avoid the risk of electrical shock, use only properly grounded equipment. Shock hazards exist if the power supply is not properly grounded. Grounding reliability can only be achieved when equipment is connected to a receptacle marked "Hospital Only" or "Hospital Grade" or the equivalent. The grounding wire must not be removed or defeated.
- To avoid the risk of electrical shock and fire hazard, inspect the power supply, plugs, cables, and power cords on a regular basis for damage.
- To avoid the risk of electrical shock, use only accessories and peripherals recommended by SonoSite, including the power supply. Connection of accessories and peripherals not recommended by SonoSite could result in electrical shock. Contact SonoSite or your local representative for a list of accessories and peripherals available from or recommended by SonoSite.
- To avoid the risk of electrical shock, use commercial grade peripherals recommended by SonoSite on battery power only. Do not connect these products to AC mains power when using the system to scan or diagnose a patient/subject. Contact SonoSite or your local representative for a list of the commercial grade peripherals available from or recommended by SonoSite.

WARNING: To avoid possible electrical shock or electromagnetic interference, verify proper operation and compliance with relevant safety standards for all equipment before clinical use. Connecting additional equipment to the ultrasound system constitutes configuring a medical system. SonoSite recommends verifying that the system, all combinations of equipment, and accessories connected to the ultrasound system comply with JACHO installation requirements and/or safety standards such as AAMI-ES1, NFPA 99 OR IEC Standard 60601-1-1 and electromagnetic compatibility standard IEC 60601-1-2 (Electromagnetic compatibility), and are certified according to IEC Standard 60950 (Information Technology Equipment (ITE)).

See the ultrasound system user guide for information and a declaration concerning Electromagnet Compatibility (EMC) for the ultrasound system.

Equipment Safety

Caution: To avoid damaging the mini-dock, do not perform high current grounding impedance testing involving this part. Accessible metal of the mini-dock is not protectively earthed. For an explanation on the grounding of the mini-dock, see SonoSite Technical Bulletin 003/06 at www.sonosite.com. Select **Support and Service** and then select your system type.

Excessive bending or twisting of cables can cause a failure or intermittent operation.

Improper cleaning or disinfecting of any part of the accessories and peripherals can cause permanent damage. For cleaning and disinfecting instructions, see [“Cleaning and Disinfecting the MDSe”](#) on page 25 and the manufacturer’s operating instructions.

Biological Safety

WARNING: Non-medical (commercial) grade peripheral monitors have not been verified or validated by SonoSite as being suitable for diagnosis.

Electromagnetic Compatibility (EMC)

See the safety section of the ultrasound system user guide and supplemental guides for information on electromagnetic compatibility with IEC 60601-1-2, including the Manufacturer’s Declaration.

MDSe Features

The MDSe provides a mobile work platform and provides a storage area for the transducers and other supplies. The MDSe includes the mini-dock with connections to accessories and peripherals. Some cables are pre-routed on the MDSe. Other cables and cords must be installed manually. See ["Cables, Connectors, and Mini-Dock Symbols"](#) on page 17.

WARNING: To avoid the risk of electrical shock or injury, the equipment shall be connected to a center-tapped single phase supply circuit when users in the United States connect the equipment to a 240V supply system.

To avoid potential injury, do not touch the squared edges on the upper columns. The squared edges may be sharp.


Figure 1 MDSe: Front and Back Views

Table 1: MDSe System Features

Number	Description
1	Ultrasound system
2	Transducer, connector, and gel storage on both sides of MDSe
3	Side panels
4	Transducer storage
5	Height adjustment pedal
6	Locking wheels
7	Locking lever
8	Black-and-white printer bay
9	VCR, DVD recorder, or color printer bay
10	Mini-dock connection panel (See “Cables, Connectors, and Mini-Dock Symbols” on page 17.)
11	Cable cover
12	Upper columns
13	Storage
14	System power cord storage knobs
15	Back panel

-
- Adjust height**
- 1 Press down on the height adjustment pedal and grasp the side panels.
 - 2 Pull up on the side panels to raise and push down to lower.
-

Lock or unlock wheels Press down on the locking levers on the front wheels to lock and push up to unlock.

Insert system	1	Turn off the ultrasound system. The mini-dock may not function properly if you connect the system while the power is on.
	2	Raise the system locking lever.
	3	Slide the system into the MDSe.
	4	Lower the locking lever to secure the system.
<i>Note:</i> If at any time the ultrasound system should lose connectivity, turn off the system and reattach mini-dock and system.		

Remove system	1	Turn off the ultrasound system.
	2	Close the lid on the ultrasound system.
	3	Lift up the locking lever and remove the ultrasound system.


Connecting Ultrasound System

Connect AC power	1	Lock the wheels and remove back panel.
	2	Ensure that the pre-routed DC power cord (A) is connected to the mini-dock and that the pre-installed AC power cord (I) is connected to the power strip.
	3	Connect the system AC power cord (K) to a hospital-grade electrical outlet.

Connecting Printer and DVD Recorder or VCR

WARNING:	To avoid the risk of electrical shock or injury, an isolation transformer must be installed on the MDSe when the following conditions exist: • More than two accessories or peripherals are used in any combination • Both color and black-and-white printers are used • Mitsubishi VCR is used with any other accessory or peripheral Contact SonoSite or your SonoSite representative for specific information on the combinations of accessories and peripherals that require installation of an isolation transformer. To avoid potential injury, do not touch the squared edges on the upper columns. The squared edges may be sharp.
Caution:	To avoid damaging the system, use only peripherals recommended by SonoSite with the ultrasound system.

The MDSe provides connections for a printer and a DVD recorder or VCR. See [Figure 2](#). For specific information on using the DVD recorders, VCR, and printers, including warnings and cautions, see the manufacturer's operating instructions.

* If Sony VCR, set dip switches: 1-4 Down and 5,6 Up.

If Panasonic VCR, set switch to Line.

** If Sony DVD recorder, use extension cable B1 and adapter B2 with B. Use adapter G2 and audio extension cables E1 with E.

If Panasonic DVD recorder, use adapter G with E.

If Panasonic VCR, use adapter F with E.

Figure 2 Connectivity of System, Printer, and a DVD Recorder or VCR

DVD Recorder and Printer

- | | |
|-----------------------|---|
| Connect cables | <ol style="list-style-type: none">1 Lock the wheels and remove back panel.2 Ensure that the pre-routed DC power cord (A) is connected to the mini-dock and that the pre-installed AC power cord (I) is connected to the power strip.3 Connect the DVD recorder cables. See Figure 2 on page 7 and “Cables, Connectors, and Mini-Dock Symbols” on page 17.<ul style="list-style-type: none">• AC power cord (I)• Pre-routed composite video cable (C)• Pre-routed S-video cable (B). For Sony, attach S-video adapter (B2) to extension cable (B1) and to S-video cable (B).• Pre-routed RS-232 cable (E):<ul style="list-style-type: none">• For Sony, attach DVD adapter (G2) to the E RS-232 connector, and attach audio extension cables (E1) to the red and white audio left/right input/output connectors.• For Panasonic, attach DVD adapter (G) to the E RS-232 connector.<i>Note:</i> RS-232 cable (E) is pre-routed through the MDSe column. You connect one end to the mini-dock. The other end has five connectors that you connect to the DVD recorder.4 Connect the printer cables. See Figure 2 on page 7 and “Cables, Connectors, and Mini-Dock Symbols” on page 17.<ul style="list-style-type: none">• Pre-routed print control cable (H)• Printer AC power cord (I)• Composite video cable (D)5 Connect system power cord (K) to a hospital-grade electrical outlet.6 Install the back panel (vents at the top) by sliding panel up and in. |
|-----------------------|---|
-

- | | |
|----------------------|--|
| Set up system | <ol style="list-style-type: none">1 Press the Setup key on the ultrasound system.2 Select Connectivity.3 In the Printer list, select desired printer.4 In the Video mode list, select the appropriate video format: NTSC or PAL.5 In the Serial Port list, select DVD.6 Restart the system to activate DVD recorder connectivity. |
|----------------------|--|

Note: Because the bar code scanner, DVD recorder, and VCR require the same port on the mini-dock, you can connect only one of these peripherals at a time.

Set up DVD recorder for recording (Panasonic)

- 1 Turn on the DVD recorder.
- 2 Press the Function key (**f**) then **V** on the system.
- 3 Press **Shift** and **Menu** on the DVD recorder.
The DVD menu appears on the system.
- 4 Select **Setup** and then press **Enter**.
A help diagram is displayed on the screen indicating how to use the DVD buttons.
- 5 From the **Remote** menu, select **Baud Rate** and **9600**.
- 6 From the **Remote** menu, select **Protocol** and **VTR1**.
- 7 From the **Picture** menu, select **Input Select** and **V (Video)**.
- 8 From the **Connection** menu, select **TV System** and **NTSC** or **PAL**.
- 9 Press **Return** on the DVD recorder to exit.
- 10 On the system, press the Function key (**f**) then **V** to exit.

Record images *Note:* Do a test recording before a patient exam to ensure that all connectivity is correct.

- 1 Turn on the DVD recorder.
- 2 Insert a DVD.
- 3 Press the **Record** key on the ultrasound system, or use controls on the DVD recorder to start and stop recording.
Each time the Record key is pressed, a new file is saved to the DVD.
(Panasonic Only) The DVD can be viewed on a PC if finalized on the DVD recorder.

Play back images

- 1 Press the Function key (**f**) then **V** on the system.
- 2 Press **Play** on the DVD recorder. Follow the manufacturer's instructions for playing back images.

VCR and Printer

Connect cables	1	Lock the wheels and remove back panel.
	2	Ensure that the pre-routed DC power cord (A) is connected to the mini-dock and that the pre-installed AC power cord (I) is connected to the power strip.
	3	Connect the VCR cables. See Figure 2 on page 7 and “ Cables, Connectors, and Mini-Dock Symbols ” on page 17. • Pre-routed S-video cable (B) • Pre-routed composite video cable (C) • Pre-routed RS-232 cable (E). For Panasonic, attach VCR adapter (F) to E. • AC power cord (I)
	4	Connect the printer cables. See Figure 2 on page 7 and “ Cables, Connectors, and Mini-Dock Symbols ” on page 17. • Pre-routed print control cable (H) • Printer AC power cord (I) • Composite video cable (D)
	5	Connect system power cord (K) to a hospital-grade electrical outlet.
	6	Install the back panel (vents at top) by sliding panel up and in.
Set up system	1	Press the Setup key on the ultrasound system.
	2	Select Connectivity .
	3	In the Printer list, select the desired printer.
	4	In the Video Mode list, select the appropriate video format: NTSC or PAL .
	5	In the Serial Port list, select VCR .
	6	Restart the system to activate VCR connectivity.
<i>Note:</i> Because the bar code scanner, DVD recorder, and VCR require the same port on the mini-dock, you can connect only one of these peripherals at a time.		
Record images	1	Ensure that the cables are correctly installed.
	2	Turn on the VCR.
	3	Ensure that the input switch is correctly set. (See the manufacturer’s instructions for location of switches.) • Sony: Dip switches: 1-4 Down and 5, 6 Up • Panasonic: Line • Mitsubishi: Dip switches: pre-set at factory
	4	Press the Record key or use the controls on the VCR to start and stop recording.
Play back images	<i>Note:</i> Do a test recording before a patient exam to ensure that all connectivity is correct.	
	1	Turn on the VCR.
	2	Press the Function key (f) then V on the system.
	3	Press Play on the VCR.

Connecting Printer Only

The MDSe provides connections for a black-and-white or color printer. See [Figure 3](#) For specific information on using the printer, including warnings and cautions, see the manufacturer's instructions.


Figure 3 Connectivity of System and a Printer on MDSe

Connect cables

- 1 Lock the wheels and remove back panel.
If installing an older black-and-white or color printer, you may need to move padded screws to back holes.
- 2 Ensure that the pre-routed DC power cord (**A**) is connected to the mini-dock and that the pre-installed AC power cord (**I**) is connected to the power strip.
- 3 Connect the cables. See [Figure 3](#) on page 11 and [“Cables, Connectors, and Mini-Dock Symbols”](#) on page 17.
 - Pre-routed print control cable (**H**)
 - Printer AC power cord (**I**)
 - Pre-routed composite video cable (**C**)
- 4 Connect system power cord (**K**) to a hospital-grade electrical outlet.
- 5 Install the back panel (vents at the top) by sliding panel up and in.

Set up system

- 1 Press the **Setup** key on the ultrasound system.
- 2 Select **Connectivity**.
- 3 In the **Printer** list, select the desired printer.
Only the printers listed are recommended for use with the system.

Print images *Note:* Do a test printing before a patient exam to ensure that all connectivity is correct.

- 1 Ensure that the cables are correctly installed.
- 2 Turn on the printer.
- 3 (Color Printer Only) Change the input setting from S-video to video.
- 4 Press the **Print** key or use the controls directly on the printer.

Connecting External Display

The MDSe provides connections for an external DVI or VGA display. For specific information on using the external display, including warnings and cautions, see the manufacturer's operating instructions.


Figure 4 Connectivity of External DVI Display on MDSe

- Connect cables (DVI display)**
 - 1 Ensure that the pre-routed DC power supply cord (A) is connected to the mini-dock and that the pre-installed AC power cord (I) is connected to the power strip.
See [Figure 4](#) and [“Cables, Connectors, and Mini-Dock Symbols”](#) on page 17.
 - 2 Remove both protective covers from the DVI cable (T).
 - 3 Connect the following cables:
 - DVI cable (T)
 - AC power cord (J)
 - DC power cord (S)
 - 4 Feed the system AC power cord (K) through the power cord protector hole at back of the lower shelf, and then connect it to a hospital-grade electrical outlet.
A firm push may be required to fit cord through protector hole.


Figure 5 Connectivity of External VGA Display on MDSe

-
- Connect cables 1 (VGA display)** Connect the cables. See [Figure 5](#) and [“Cables, Connectors, and Mini-Dock Symbols”](#) on page 17.
- External display VGA cable (**M**) to the external display and mini-dock
 - External display power cord (**L**)
- 2 Connect the system power cord (**K**) to a hospital-grade electrical outlet.

Turn on display and adjust angle

Adjust the display to the proper angle prior to use.

- 1 Turn on the display. (The switch is on the back.)
- 2 Turn on the ultrasound system.
- 3 Holding the display with both hands, rotate it forward or backward to the desired angle. Additional pressure may be required.
- 4 (VGA only) On the display, press + or – to select the desired video mode.


Figure 6 Adjust the Angle of the External Display

Caution:

A set screw controls the tightness of the display angle. To avoid damage to the display, do not adjust the set screw. Turning the set screw loosens the display, causing the display to drop.

Connecting Triple Transducer Connect

The Triple Transducer Connect is a three-transducer module for the MDSe that allows connection of up to three transducers simultaneously while the ultrasound system is docked.

- | | |
|---|---|
| Connect transducer connector to system | <ol style="list-style-type: none">1 Raise locking lever.2 Pull the transducer latch up and rotate 90°.3 Align and insert the transducer connector to the connector on the bottom of the ultrasound system.4 Turn the latch back 90° and press in, securing the transducer connector to the ultrasound system.5 Place the system into the MDSe. <p>Connection is made when locking lever is pushed down.</p> |
|---|---|


Figure 7 Connect Transducer Connector to System

- | | |
|--|--|
| Connect transducer to Triple Transducer Connect | <ol style="list-style-type: none">1 Pull the transducer latch up and rotate 90°.2 Align and insert the transducer connector to the connector on the bottom of the Triple Transducer Connect.3 Turn the latch and press down, securing the transducer connector to the Triple Transducer Connect. |
|--|--|

Select transducer *Note:* When you use the Triple Transducer Connect, there may be a minor reduction of transducer penetration.

Push the transducer select button on the Triple Transducer Connect to select the desired transducer. See [Figure 7](#) on page 16.

- When selecting another transducer, the system restarts. The light on the button flashes momentarily and then turns to a solid green.
- Only one transducer is accessible at a time.

Cables, Connectors, and Mini-Dock Symbols

The cable kit includes cables for both the MDSe and the MDS Lite. You need to use only the cables for your product.

Table 2: Connectivity Symbols on Mini-Dock

Symbol	Definition
	DC input
	Print control
	USB
	Ethernet
	RS-232 (DVD recorder, VCR, or bar code scanner)
	S-video out
	S-video in
	DVI video out
	Composite video out
	Audio out
ECG	ECG/Footswitch

Table 3: Pre-routed Cables









Letter	Name	Illustration	Connector 1 [Illustration] [Connects to]	Connector 2 [Illustration] [Connects to]
A	DC power supply cord 6.8 ft./2 m		 Mini-dock	 Power supply
B	S-video cable (playback) 6 ft./1.8 m		 Mini-dock	 DVD recorder or VCR
C	Composite video cable 6 ft./1.8 m		 Mini-dock	 DVD recorder or VCR
E	RS-232 cable (control/audio) 6.5 ft./2 m		 Mini-dock	  DVD recorder or VCR
H	Print control cable 6 ft./1.8 m		 Mini-dock	 Printer

Table 4: Manually Installed Cables and Connectors

Letter	Name	Illustration	Connector 1 [Illustration] [Connects to]	Connector 2 [Illustration] [Connects to]	Product Used On
B	S-video cable (playback) 6 ft./1.8 m		 Mini-dock	 DVD recorder or VCR	MDS Lite
B	S-video cable (record) 6 ft./1.8 m		 Mini-dock	 DVD recorder or VCR	MDS Lite
B1	S-video extension cable (Sony DVD recorder only) 12 in./30.5 cm		 S-video adapter	 DVD recorder	MDSe
B2	S-video adapter (Sony DVD recorder only)		 S-video cable	 S-video extension cable	MDSe
C	Composite video cable 6 ft./1.8 m		 Mini-dock	 DVD recorder or VCR	MDS Lite
D	Composite video cable 12 in./30.5 cm		 Printer	 DVD recorder or VCR	MDSe

Table 4: Manually Installed Cables and Connectors










Letter	Name	Illustration	Connector 1 [[Illustration] [Connects to]	Connector 2 [[Illustration] [Connects to]	Product Used On
E	RS-232 cable (control/audio) 6.5 ft./2 m		 Mini-dock	 DVD recorder or VCR	MDS Lite
E1	Audio extension cables (Sony DVD recorder only) 12 in./30.5 cm		 RS-232 cable (audio output/input)	 DVD recorder	MDSe
E1	Audio extension cables (Sony DVD recorder only) 12 in./30.5 cm		 RS-232 cable (audio output/input)	 DVD recorder	MDSe
F	VCR adapter (Panasonic only) 12 in./30.5 cm		 RS-232 cable	 VCR	MDSe MDS Lite
G	DVD adapter (Panasonic only) 12 in./30.5 cm		 RS-232 cable	 DVD recorder	MDSe MDS Lite

Table 4: Manually Installed Cables and Connectors


















Letter	Name	Illustration	Connector 1 [Illustration] [Connects to]	Connector 2 [Illustration] [Connects to]	Product Used On
G2	DVD adapter (Sony only) 12 in./30.5 cm		 RS-232 cable	 DVD recorder	MDSe MDS Lite
H	Print control cable 6 ft./1.8 m		 Mini-dock	 Printer	MDS Lite
I	AC power cord 18 in./45.7 cm		 Power supply	 Power strip	MDSe MDS Lite
I	AC power cord 18 in./45.7 cm		 DVD recorder or VCR	 Power strip	MDSe MDS Lite
I	AC power cord 18 in./45.7 cm		 Printer	 Power strip	MDSe MDS Lite
J	AC power cord 39 in./1 m		 Printer	 Power strip	MDS Lite
K	System AC power cord 10 ft./3 m		 Power strip	 Wall outlet	MDSe MDS Lite (ships in separate box)
L	AC power cord 6 ft./1.8 m		 External display	 Power strip	MDSe

Table 4: Manually Installed Cables and Connectors

Letter	Name	Illustration	Connector 1 [[Illustration] [Connects to]	Connector 2 [[Illustration] [Connects to]	Product Used On
M	External display VGA 3 ft./ .9 m		 Mini-dock	 Display	MDSe
S	External display power supply		 Display	 Power supply	MDSe
T	External display DVI 3 ft./ .9 m		 Mini-dock	 Display	MDSe

Troubleshooting

Symptom	Solution
External VGA display does not show an image.	• Check the video format setting in system setup. See “Turn on display and adjust angle” on page 15. • Check the cable connections on display and stand. See Figure 4, “Connectivity of External DVI Display on MDSe” on page 13 and “Connectivity of System, Printer, and a DVD Recorder or VCR” on page 7. • Check the display to ensure that it is turned on and set up properly (i.e., ready to show correct video format: S-video, composite, or RGB). See the display manufacturer’s instructions, if necessary.
External display does not project sound.	• Check the cable connections on display, stand, and VCR. See Figure 4, “Connectivity of External DVI Display on MDSe” on page 13 and “Connectivity of System, Printer, and a DVD Recorder or VCR” on page 7.
DVD recorder does not record.	• Check the serial port setting in system setup, including Record key setting and DVD recorder setup. See “Set up system” on page 8 and “Set up DVD recorder for recording (Panasonic)” on page 9. • Check the cable connections on stand and DVD recorder. See Figure 2, “Connectivity of System, Printer, and a DVD Recorder or VCR” on page 7. • Check the DVD recorder to ensure that it is turned on and set up properly. See the manufacturer’s instructions, if necessary.
DVD recorder does not play back.	• Check the cable connections on display, stand, and DVD recorder. See “Connectivity of External DVI Display on MDSe” on page 13 and Figure 2, “Connectivity of System, Printer, and a DVD Recorder or VCR” on page 7. • Check the display setting to ensure that it is set to S-video. Press the +/- key to toggle video format on the display (default setting is Composite).

Symptom	Solution
VCR does not record.	• Check the video format setting in system setup, including Print/VCR key setting. See “Set up system” on page 10. • Check the cable connections on stand and VCR. See Figure 2, “Connectivity of System, Printer, and a DVD Recorder or VCR” on page 7. • Check the VCR to ensure that it is turned on and set up properly. See the VCR manufacturers’ instructions, if necessary. • Check the input switch setting. See “Record images” on page 10.
VCR does not play back.	• Check the cable connections on display, stand, and VCR. See “Connectivity of External DVI Display on MDSe” on page 13 and Figure 2, “Connectivity of System, Printer, and a DVD Recorder or VCR” on page 7. • Check the display setting to ensure that it is set to S-Video. Press the +/- key to toggle video format on the display (default setting is Composite). • Check the input switch setting. See “Record images” on page 10.
Printer does not work.	• Check the printer selection in system setup. See “Set up system” on page 11. • Check the cable connections. See “Connecting Printer and DVD Recorder or VCR” on page 6 or “Connecting Printer Only” on page 11. • Check the printer to ensure that it is turned on and set up properly. See the printer manufacturer’s instructions, if necessary.

Cleaning and Disinfecting the MDSe

WARNING:	To avoid electrical shock, remove the system and extra battery from the MDSe and unplug from power supply before cleaning.
	To avoid injury, always use protective eyewear and gloves when performing cleaning and disinfecting procedures.
	To avoid infection, if a pre-mixed disinfectant solution is used, observe the solution expiration date, and ensure that the date has not passed.
	To avoid infection, the level of disinfection required for a product is dictated by the type of tissue it contacts during use. Ensure the solution strength and duration of contact are appropriate for the equipment. For information, see the disinfectant label instructions and the recommendations of the Association for Professionals in Infection Control and Epidemiology (APIC) and FDA.
Caution:	Use only recommended cleaners or disinfectants on surfaces. Immersion-type disinfectants are not tested for use on surfaces.

Clean and disinfect	The exterior surface of the MDSe can be cleaned and disinfected using a recommended cleaner or disinfectant.
	1 Remove the system, unplug the power supply, and detach any cables.
	2 Clean the surfaces using a soft cloth lightly dampened in a mild soap or detergent cleaning solution.
	Apply the solution to the cloth rather than the surface.
	3 Wipe the surfaces with the disinfecting solution. Theracide disinfectant is recommended.
	4 Air dry or towel dry with a clean cloth.

Cleaning the Peripherals

For information about cleaning the peripherals, see the manufacturer's operating instructions.

Specifications

MDSe Dimensions

Width: 22.6 in. (57.40 cm)
Length: 31.9 in. (81.03 cm)
Height: 34.2 - 46.2 in. (86.88–117.35 cm)
Height range: 12 in. (30.48 cm)
Weight: 78 lbs. (39.92 kg)

Temperature and Humidity Limits

Operating: 10–40°C; 15–95% R.H.
Shipping/Storage: -35–65°C; 15–95% R.H.

Peripherals

See the manufacturer's specifications for the peripherals.

Medical Grade

Black-and-white printer

Recommended sources for printer paper: Contact Sony at **1-800-686-7669** or **www.sony.com/professional** to order supplies or to obtain the name of the local distributor.

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